

# Instruction Manual for Microprocessor Centrifuge DSC-302SMD (DC brushless motor)





### Warnings:

- \* Do not open the lid before the centrifuge stops completely.
- \* Do not remove the centrifuge during operation.
- \* Do not put objects on the top of the lid.
- \* Do not disassemble the centrifuge by yourself. In the event of a malfunction, the centrifuge should be repaired by a qualified technician. Before repairing, turn off the power.

#### Steps:

- 1. Place the centrifuge on a flat surface. Plug the centrifuge into a properly grounded outlet. Turn on the main power switch on the right side of the centrifuge.
- 2. Push the button , and then the screen is illuminated. Press the button to open the lid. The icon is flash on the screen. Insert an even number of tubes in a symmetrical position into the rotor.
- 3. Put down and press down the lid. Make sure the lid is closed well. The icon will disappear on the screen.
- 4. Push the button pred, the speed value is flash one the screen. Push the button or to set the speed value. Push the button or any other function button to enter the speed value.

- 5. If you want to set RCF value instead of setting Speed value, please push the button the button the button to set RCF value instead of setting Speed value, please push the button to speed and then push the button to set RCF value. Push the button to set RCF value. Push the button to set RCF value.
- 6. Push the button , the time value is flash on the screen. Push the button or to set seconds. Push again to set minutes. Push the button or any other function button to enter the time value.
- 7. Push the button ROTOR. The rotor type value is flash. Push or to set rotor type (refer to the below table.) Push the button ROTOR again. The tube adapter size value is flash. Push the button or to choose the tube adapter size. Push the button ROTOR to enter the rotor value.

Rotor type	Rotor type on the screen	Tube adapter size
AR-1512	Ar 12	2ml, 7ml, 15ml
AR-1524	8-24	2ml, 7ml, 15ml
UAR-1512	UA- 12	2ml, 7ml, 15ml
AR-5006	Rr 05	50ml
SB071550	-56-	7ml, 15ml, 50ml

- 8. After setting parameters, push the button to start the centrifuge. During operation, push the button to stop the centrifuge.
- 9. Push the button start the centrifuge directly. Release the button and the centrifuge will stop immediately. The centrifuge will run to the current speed setting value.

- 10. Brake setting: the initial setting is auto brake and the screen will show **BRK**. If you want to cancel the auto brake function (which means when the machine stop rotating, there is no auto brake function, and the machine will stop rotating naturally.) Push the button continuously, and at the same time, push the button the BRK will disappear. Repeat the same step, the auto brake function can restore.
- 11. Using memory program setting:
  - 11.1 How to use memory program: push the button , the icon is flash on the screen. Push the button or to choose the memory number 1,2,3, ..., 9. Then, push the button start the centrifuge
  - 11.2 Parameter setting: Push the button or to choose the memory number 1,2,3, ..., 9. Push the button to set the speed value.

    Push the button to set the time value. Push the button to set rot and tube adapter size.

    Push the wew button to finish the memory setting.
  - 10.3 How to exit the memory mode: Push the button . Push the button or , until there is no number beside the icon . (Note: Another way to exit the memory mode. Push the button to turn off the power, and then push the button again to exit the memory mode.)
- 12. Under the conditions of power failure or the open function failure, please insert the "Emergency Bar" into the hole at right hand side of centrifuge to force the lid open.

Other functions: (Under no operation condition, it can execute other functions.)


- 1. Adjust the contrast of the characters on the screen:
  - 1.1 Push the button , and do not release. At the same time push the button to increase the contrast. The icon ACC is flash on the screen.
  - 1.2 Push the button , and do not release. At the same time push the button to decrease the contrast. The icon DEC is flash on the screen.
- 2. Adjust the brightness of the background of the screen:
  - 2.1 Push the button , and do not release. At the same time push the button brightness. The icon ACC is flash on the screen.
  - 2.2 Push the button , and do not release. At the same time push the button to decrease the

brightness. The icon DEC is flash on the screen.

### Maintenance and care:

- Keep the centrifuge clean all the time. Cleaning can be done with a damp cloth. Avoid the use of solvents as they may damage the product housing, rotor.
- 2. Check if the rotor is in the center. In case, the rotor is slanted. Stop using machine and contact the distributor for repair.

### **Circuit Diagram:**



## **Service Manual**

# **Troubleshooting Guide**

Caution: Always disconnect the power cord before troubleshooting.

Trouble / error message		age	Cause	Remedy
Motor is not working	X	<b>∌</b> £₽		Turn the main power switch off, and then turn it on again.
			1. Rotor type setting is not correct.	Please check and re-set the rotor.
			2. Foreign bodies block the motor.	2. Sweep away foreign bodies.
			3. Motor is defective.	3. Contract the distributor for repair.
Over heat	X	<b>● EF</b>		Turn the main power switch off, and
	•	2.5		then turn it on again.
			1. Rotor type setting is not correct.	1. Please check and re-set the rotor.
			2. Over time use.	2. Stop using to let the centrifuge cool
				down. And then restart it.
			3. Motor is defective.	3. Contract the distributor for repair.
Over current	X	EA		Turn the main power switch off, and
				then turn it on again.
			1. Rotor type setting is not correct.	1. Please check and re-set the rotor.
			2. Rotor is overloaded $\circ$	2. Reduce the loading on the rotor.
			3. Motor is defective.	3. Contract the distributor for repair.
Over voltage	X	® Ei		Turn the main power switch off, and
				then turn it on again.
			1. Rotor type setting is not correct.	1. Please check and re-set the rotor.
			2. Switching power is defective.	2. Contract the distributor for repair.
	X	EL	1. When the motor is rotating, push	1. Waiting for motor stop completely,
			the button 😃 .	and then push the button OPEN,
				and then close the lid to restart the
				centrifuge.
Unbalance	<b>4</b>	æ E∐	1. Rotor is unbalanced.	1. Open the lid and set the rotor to be
				balanced.
			2. The surface which the centrifuge is	2. Place the centrifuge on a flat
_			on is not flat.	surface.
Lid open error		> Ed	1. During operation, the lid is open.	1. Waiting for the motor stop
				completely. And close the lid.
Speed error	<b>飛</b> (	<b>®</b> Er	1. Rotor type setting is not correct.	Open the lid. Check and reset the rotor.
			2. Motor, PC board, or switching	2. Contract the distributor for repair.
			power is defective, so the speed	
			can't reach to the setting value.	
No display on the scree	en		1. Power is off.	1. Turn on the power.

	2. Dead power output $\circ$	2. Change to a different output.
	3. Burnt fuse	3. Replace a new fuse.
	4. PC board is defective.	4. Contract the distributor for repair.
Burnt fuse	1. Fuse is wasted	1. Replace a new fuse.
	2. Switching power is defective.	2. Contract the distributor for repair.
Lid can't be open	1. Lid lock is defective.	1. Contract the distributor for repair.
	2. PC board is defective.	2. Contract the distributor for repair.
	3. PC board is crashed.	3. Turn off the main power switch, and
		turn on the power again.
Noise is too loud.	1. Rotor is unbalanced.	1. Stop the machine. Open the lid
		and set the rotor to be balanced.
	2. Rotor type setting is not correct.	2. Please check and re-set the rotor.

Note: To clear the error message **EP**, **EE**, **EP**, **EU**, turn the main power switch off, and turn on it again.

For the other error messages, push the button open to clear them.