

Accessories and Light source for **STEREO MICROSCOPES**


- * The cold light source has high intensity fiber optic illumination and accessories that allow you to position the light for the best possible viewing.
- * The 150W lamp is equipped with a variable control.
- * Both power suppliers are designed for easy lamp replacement.



Type A power supplier
Power: 110V/220V switchable
Size: 280x120x170mm
Housing material: metal
Equipped with automatic power off device in lamp housing.



Type B power supplier
Power: 85V~265V switching power
Size: 265x170x135mm
Housing material: ABS plastic
Adopt over loading protect set on PC board.



Fibers and accessory:
(a) Koehler illumination attachment **KA-9** : Attached on Single or Dual light guide to adjust and obtain even illumination.
(b) Single light guide **SLG-500** : Rigid flexional tube, length 500mm
(c) Dual light guide **DLG-500** : Rigid flexional tube, length 500mm
(d) Ring light guide **RLG-1000** : Flexible tube, length 1000mm

Easy order information:

	Type A Power Supplier	Type B Power Supplier
Single Light Guide SLG-500	SLG-500A	SLG-500B
Dual Light Guide DLG-500	DLG-500A	DLG-500B
Ring Light Guide RLG-1000	RLG-1000A	RLG-1000B

Note: KA-9 is optional for Model SLG-500A, DLG-500A, SLG-500B and DLG-500B.



LED-60
Adjustable LED ring illuminator



FRL-8
8W fluorescent ring illuminator



MS-H26
Stage size: 180x155x26mm
Moving range: 75x55mm with graduation of 0.1mm
Connecting size: Diam. 95mm



MS-H43
Stage size: 180x155x43mm
Moving range: 75x55mm with graduation of 0.1mm
Connecting size: Diam. 95mm
Glass plate: Diam. 95mm & 360 degree rotatable



HL-20
Adjustable 12V20W halogen lamp