CENTRIFUGE



TABLE TOP CENTRIFUGE

DSC-200T & DSC-200A series







DSC-200T

DSC-200A-1

DSC-200A-2

Features:

- Small centrifuge can save your working desk space.
- ASB outer shell can resist collision.
- With safety switch, when lid is lifted during operation, the rotor will stop spinning.
- Quality certificate: ISO 9001, ISO 13485.

- There are three models for your choices.
- Auto balancing system can offer quiet and stable operation.
- All models can accommodate variable centrifuge tube size from 1.5ml to 15ml without changing rotor.
- Model DSC-200A-2 has a tachometer which shows both speed and R.C.F value.

Specifications:

Model	Max. Speed (rpm)	Timer (Min.)	Speed & R.C.F. display	Manual brake	Motor	Weight (kg) (N.W./G.W.)	Dimension (mm)(WxDxH)
DSC-200T	3,000 at 60Hz 2,500 at 50Hz	0~30 min.	None	None	Brushless induction	3.5 / 4.5	280x280x255
DSC-200A-1	6,000	mechanical timer	None	Yes	DC	4.5 / 5.5	280x280x265
DSC-200A-2	6,000	unioi	Yes	Yes	DC	4.5 / 5.5	280x280x265
Dower avealy A C440V 60U = (220V F0/60U = aveilable)							

Power supply: AC110V 60Hz (220V 50/60Hz available).

Specifications of rotors:

Model	Speed (rpm)	Туре	Max. R.C.F. (g)	Capacity	Standard accessory		
Rotor for centrifuge DSC-200T							
AR-1506	300~3.000	45° angle	1,066	6x(10~15)ml	(10~15)ml tube adapter		
AR-1508	300~3,000		1,066	8x(10~15)ml			
Rotor for centrifuge DSC-200A-1 and DSC-200A-2							
UAR-1508	300~6,000	45° angle	4,265	8x(10~15ml)	(10~15)ml tube adapter		
Note: Tube adapters for (1.5~2.0)ml or (5~7)ml tubes are optional.							

Easy order information:

Order No.	Item
DSC-156T	Centrifuge DSC-200T with rotor AR-1506
DSC-158T	Centrifuge DSC-200T with rotor AR-1508
DSC-N158T	Centrifuge DSC-200A-1 with rotor UAR-1508
DSC-N158A	Centrifuge DSC-200A-2 with rotor UAR-1508
CTA15	Tube adapter for (10~15)ml tube; inner size: Φ17.5 x L98mm
CTA7	Tube adapter for (5~7)ml tube; inner size: Φ13.7 x L69.5mm
CTA2	Tube adapter for (1.5~2.0)ml tube; inner size: Φ11.3 x L34mm





CTA7 CTA15

No.1, Lane 393, Fu De 1 Road, Xizhi Dist., New Taipei City 22149, Taiwan T:886-2-26942339.26945547.26941738 F:886-2-26941695 ⊠ digilab @ ms23.hinet.net www.digisystemlab.com