

## ST-600 series

# STEREO MICROSCOPES

### Feature:

- \* The objective is 2X/4X of turret style.
- \* The image is sharp and clear over a wide field of view.
- \* To avoid killing the cell due to high temperature of lamp, a fluorescent lamp for 10,000 hours lighting is available.
- \* The microscope frame is made of aluminum alloy with stable structure.
- \* It is excellent optics for school teaching and industry manufacturing.



**ST-604**



**ST-606 with digital camera**

### Specifications:

Model	ST-604	ST-606
Observation tube	Binocular head inclined at 45 degree	Trinocular head inclined at 45 degree
Focusing	Rack and pinion focusing assembly	
Stage	220mm × 155mm	
Eyepiece	WF10X	
Objective	2X / 4X of turret style	
Working Distance	90mm	
Total Magnification	20X / 40X (The combination of different eyepieces and supplemental lenses yield magnification as high as 160X)	
Illumination	Incident: 12V20W halogen lamp, Transmitted: 12V10W halogen lamp or 5W fluorescent lamp	
Power Supply	110V (220V available)	
Standard Accessories	Dust cover, paired eyecaps, frosted glass insert, white/black insert	
Optional Accessories	WF15X or WF20X eyepiece, 0.5X, 0.75X, 1.5X or 2X supplemental lenses, 150W cold light source, LED ring illuminator, 8W fluorescent ring illuminator, hi-resolution color CCD, CCD adapter, digital camera adapter	