

DIFFERENTIAL BLOOD CELL COUNTER

DBC-2 / DBC-6 / DBC-9



Features:

- * ABS outer shell can resist collision.
- * Color charts of white blood cells are supplied. (DBC-6 including BASO., EOSIN., SEG., LYMPH., MONO.; DBC-9 including BASO., EOSIN., MYELO., JUVEN., STAB., SEG., LYMPH., MONO.)
- * At one hundred increment of the total record, the bell sound "DON" will be heard (for DBC-6, DBC-9)
- * Reset knob can easily clear all figures.
- * Quality certification: ISO 9001, ISO 13485

Specifications:

Model	Total of keys	Total of windows	Totalizer	Figure range	Overall dimension (mm) (WxDxH)
DBC-2	2	2	No	0-999	100x80x50
DBC-6	5	6	Yes	0-999	250x80x50
DBC-9	8	9	Yes	0-999	320x80x50

HAND-HELD COUNTER

DBC-1



Features:

- * Metal housing.
- * Easy held by one hand.
- * Reset knob.
- * Quality certification: ISO 9001, ISO 13485

Specifications:

Model	Total of keys	Total of windows	Figure range
DBC-1	1	1	0-9999

BOILING STERILIZER

BS-1836



Features:

- * All hydraulically pressed 18.8 stainless steel, can withstand electrolysis without corrosion.
- * Long life immersion stainless nickel pipe heater, explosion & corrosion proof.
- * Precision mechanical timer with bell.
- * Waterless overheating prevention system with auto power cut-off.
- * Quality certification: ISO 9001, ISO 13485

Specifications:

Model	Timer	Heater	Thermostat	Tray dimension (mm) (WxDxH)	Chamber dimension (mm) (WxDxH)	Overall dimension (mm) (WxDxH)
BS-1836	0-60 min. mechanical timer	1500W	Temp. over to 150 degree C auto cut-off power	350x165x30	360x180x125	525x190x170