

PLATELET SHAKER

PS-48 & PS-96



PS-48



PS-96

Features:

- * Platelet shaker PS-48 and PS-96 are widely used in blood banks.
- * A stationary base drives the storage frame side-to-side 38mm.
- * Removable open mesh shelves provide air circulation.
- * A ball bearing drive motor and internal fan produce durable and reliable operation.
- * Stainless steel shelves can be easily cleaned.
- * PS-48 and PS-96 are easy to use.
- * Quality certification: ISO 9001, ISO 13485


Specifications:

Model	PS-48	PS-96
Capacity	42 bags	84 bags
Speed	70rpm (fixed)	
Number of shelf	7 shelves	
Overall dimension (mm)(WxDXH)	460x350x360	840x350x360

Power supply: AC 110V 60Hz (220V 50/60Hz available)

PLATELET ROTATOR

PR-Series



PR-6

Features:

- * Platelet rotator PR-2, PR-4, PR-6, PR-8 are widely used in blood banks.
- * Robust, practical, and durable design for your economic choice.
- * The storage frame rotates 360° smoothly and steadily.
- * Platelet bags can be loaded easily on the stainless steel basket.
- * Attached with a tachometer, users can monitor speed easily.
- * Quality certification: ISO 9001, ISO 13485

Specifications:

Model	PR-2	PR-4	PR-6	PR-8
Number of basket	1	2	3	4
Capacity	12 bags	24 bags	36 bags	48 bags
Speed	0~9 rpm			
Overall dimension (mm)(WxDXH)	370x255x290	510x255x290	650x255x290	790x255x290

Power supply: AC 110V 60Hz (220V 50/60Hz available)



DIGISYSTEM LABORATORY INSTRUMENTS INC.

NO. 1, LANE 393, FU DE 1 ROAD, XIZHI DIST., NEW TAIPEI CITY, TAIWAN
 TEL: 886-2-26941738 . 26942339 . 26945547 FAX: 886-2-26941695
 E-MAIL: digilab@ms23.hinet.net INTERNET: www.digisystemlab.com
 We reserve the rights to change the specifications without notice.